

# About Omics Group

[OMICS Group](#) International through its Open Access Initiative is committed to make genuine and reliable contributions to the scientific community. [OMICS Group](#) hosts over 400 leading-edge peer reviewed Open Access Journals and organize over 300 International Conferences annually all over the world. OMICS Publishing Group journals have over 3 million readers and the fame and success of the same can be attributed to the strong editorial board which contains over 30000 eminent personalities that ensure a rapid, quality and quick review process.

# About Omics Group conferences


- [OMICS Group](#) signed an agreement with more than 1000 International Societies to make healthcare information Open Access. [OMICS Group](#) Conferences make the perfect platform for global networking as it brings together renowned speakers and scientists across the globe to a most exciting and memorable scientific event filled with much enlightening interactive sessions, world class exhibitions and poster presentations
- Omics group has organised 500 conferences, workshops and national symposium across the major cities including San Francisco, Omaha, Orlando, Raleigh, Santa Clara, Chicago, Philadelphia, United Kingdom, Baltimore, San Antonio, Dubai, Hyderabad, Bangaluru and Mumbai.

**Lasers, Optics and Photonics 2014 – Philadelphia, USA, Sept. 9, 2014**

# **Linear optical methods as ultrashort pulse diagnostics**



**A. Börzsönyi, A.P. Kovács, K. Osvay**

*Department of Optics and Quantum Electronics, University of Szeged, Hungary  
Extreme Light Infrastructure – Attosecond Light Pulse Source, Hungary  
CE Optics Research and Development Ltd., Hungary*



# Ultrashort laser pulses

Only a few optical cycles under the pulse envelope



$$E(t, z) = \frac{1}{\sqrt{2\pi}} \int_{-\infty}^{\infty} \tilde{E}(\omega) \exp[i\omega(t - n(\omega) \cdot z/c)] d\omega.$$


Where the spectral phase is

$$\varphi(\omega, z) = \omega \cdot n(\omega) \cdot z/c$$

Fourier analysis suggests that short pulses have broadband spectra

# Phase relations of ultrashort pulses

**Carrier-envelope phase:** relation between the carrier wave and the pulse envelope.



## Use of CEP stabilized pulses


- Attosecond physics and high-harmonic generation
- High precision optical frequency and time measurements (by stable "frequency-comb")
- High precision refractive index measurements
- Calibration of astronomical mirrors

# Time-domain pulse distortions

**Dispersion of the spectral phase related strongly to the pulse shape:**

$$\varphi(\omega, z) = \varphi(\omega_0, z) + \left. \frac{d\varphi(\omega, z)}{d\omega} \right|_{\omega_0} \cdot (\omega - \omega_0) + \frac{1}{2} \left. \frac{d^2\varphi(\omega, z)}{d\omega^2} \right|_{\omega_0} \cdot (\omega - \omega_0)^2 + \frac{1}{6} \left. \frac{d^3\varphi(\omega, z)}{d\omega^3} \right|_{\omega_0} \cdot (\omega - \omega_0)^3 + \dots$$


Case of propagation in vacuum



# Space-domain pulse distortions

**Spectral components of the pulse** propagate to different directions

- Various spatiotemporal pulse deformations
- Also possible to tailor the pulse length in a controlled way: stretcher-compressor systems of chirped pulse amplifiers




# Space-domain pulse distortions


## Two different definitions of angular dispersion

- Propagation direction of spectral components
- Spectral phase fronts of spectral components
- Same for plane waves
- Different for Gaussian beams, phase front angular dispersion slowly diminishes during propagation

(a) Plane waves



(b) Gaussian beams






# Self-referenced pulse measurements

Based on frequency conversion of the autocorrelation signal


- Typical methods: FROG, SPIDER, WIZZLER, etc.
- Requires nonlinear process
- Usually difficult and iterative algorithms



# Linear optical methods

## Spectrally (and spatially) resolved interferometry




- Based on the combination of a two-beam interferometer and an imaging spectrograph
- No need for nonlinear processes
- Easy and straightforward algorithm
- High precision easily achievable
- In some cases, it is restricted to measure relative values only



# Fourier-transformation algorithm

## Evaluation is based on Fourier-domain filtering


- Record the spectrally resolved interferogram
- Keep only the relevant part in the Fourier domain
- Apply inverse Fourier-transformation
- Complex angle of the received spectrum gives the phase
- Phase derivatives calculated from polynomial fitting



# Stationary Phase Point Method

**Effective when the sample pulse is much longer than the reference**


- Can be used to characterize stretchers/compressors
- Short pulse is only coherent in a small spectral region
- Stationary phase point is where the modulation minimal
- SPP can be scanned through the spectrum by changing the delay of the reference pulse
- Delayed between the edges of the spectrum tells the length of the stretched pulse



# Spectrally and spatially resolved interferometry

## A bit more general description of the SRI


- Based on the combination of a two-beam interferometer and an two-dimensional imaging spectrograph
- A small angle between the incoming beams at the spectrograph's slit helps the algorithm and visualizes the spectral phase



# Spectrally and spatially resolved interferometry

## A bit more general description of the SRI

- Based on the combination of a two-beam interferometer and an two-dimensional imaging spectrograph
- A small angle between the incoming beams at the spectrograph's slit helps the algorithm and visualizes the spectral phase




# 2D Fourier-transformation algorithm

The 2D implementation of the Fourier-domain filtering works well

The shape of the spectral phase translates directly to the shape of the fringes


- (a) First order dispersion only (**group delay**)
- (b) Second order dispersion only (**group delay dispersion**)
- (c) Third order dispersion only (**third order dispersion**)



# Angular dispersion measurement

## Slight modification of the general SSRI setup

- The interfering beams are mirrored (e.g. one more reflection)
- Spectral density of the fringes correspond to the angle between the interfering beams
- If the beam has angular dispersion, the mirroring effect doubles the spectral dependence of the phase front direction
- The absolute phase front angular dispersion can be measured






# Alternative angular dispersion measurement

**Requires a 2D spectrograph, a beam rotation stage and a focusing element only**




- Focus the beam onto the slit of the spectrograph
- The angle of the spectral components translated to position
- Spectral components with different position means angular dispersion
- Spectrograph resolves the beam spectrally
- The tilt of the spectrogram tells the angular dispersion
- Polarization (beam orientation) rotation needed for full characterization



# Alternative 2D angular dispersion measurement


## Neither spectrograph, nor beam rotation needed

- Calibrated spectral filtering (e.g. Fabry-Perot etalon) creates well defined, separated spectral peaks
- Focusing element translates angle into position
- Spot distances gives the angular dispersion in 2D



# Phase relations of ultrashort pulses

**Carrier-envelope phase:** relation between the carrier wave and the pulse envelope.




## Use of CEP stabilized pulses

- Attosecond physics and high-harmonic generation
- High precision optical frequency and time measurements (by stable "frequency-comb")
- High precision refractive index measurements
- Calibration of astronomical mirrors

# Methods to measure CEP drift

## The f-to-2f scheme – the nonlinear way



$$f_1 = f_{CEO} + n \cdot f_{rep}$$

$$f_2 = 2 \cdot (f_{CEO} + m \cdot f_{rep})$$

$$f_{beat} = f_2 - f_1 = f_{CEO}$$


$$\varphi_{CEO} = 2\pi \cdot f_{CEO} / f_{rep}$$

### Requirements:

- Octave-broad bandwidth
- 2<sup>nd</sup> harmonic generation

# Methods to measure CEP drift


## Multiple-beam interferometer – the linear way



K. Osvay et al., Opt.Lett. **32**, 3095 (2007).

### Evaluation steps:


- (1) Record interference pattern
- (2) FFT and filter the spectrum
- (3) Inverse FFT
- (4) Complex angle gives the spectral phase difference
- (5) Fitting Taylor-series
- (6) Calculate  $CEP = \phi_0 - GD \cdot \omega_0$



L. Lepetit et al., JOSA B **12**, 2467 (1995).

# Methods to measure CEP drift



## Multiple-beam interferometer – the linear way





- **Linear**
- **Scalable**
- **Bandwidth independent**
- **UV and far infrared lasers**
- **Applicable to a wide range of lasers:**
- **(sub-)picosecond lasers**

# Cross-calibration with f-to-2f method

## Sinusoidal intracavity modification of CEP drift




## Random intracavity modification of CEP drift



# CEP drift noise in a Ti:S amplifier

## Spectrally resolved interferometry – relative CEP change





- General setup of the measurement (Mach-Zehnder interferometer)
- Method: Spectrally Resolved Interferometry (SRI)
- Main source of CEP change: Ti:S amplifier crystal
- Measuring the change in CEP between non- and amplified pulses




# CEO noise and drift measurement

- Effect of crystal cooling on CEP stability:




- Coolant temperature:



Temperature values are related to the cooling water

CEP noise:



# Further developments

- **Spectrally resolved interferometry is an efficient linear optical method of ultrashort pulse measurement**
- **Most effective if used together with self-referenced characterization techniques**

**Thank you for your attention!**

#### **Acknowledgements:**

This research was supported by the European Union and the State of Hungary, co-financed by the European Social Fund in the frameworks of TÁMOP 4.2.4. A/2-11-1-2012-0001 'National Excellence Program' and TÁMOP 4.2.2.A-11/1/KONV-2012-0060 – „Pulse lasers for use in materials science and biophotonics” and „TÁMOP-4.2.2.C-11/1/KONV-2012-0010 – „Supercomputer, the national virtual laboratory”.

**SZÉCHENYI** 2020



**Európai Unió**  
Európai Szociális  
Alap



**BEFEKTETÉS A JÖVŐBE**

# Let Us Meet Again

We welcome all to our future group conferences  
of Omics group international

Please visit:

[www.omicsgroup.com](http://www.omicsgroup.com)

[www.Conferenceseries.com](http://www.Conferenceseries.com)

<http://optics.conferenceseries.com/>