

About Omics Group

[OMICS Group](#) International through its Open Access Initiative is committed to make genuine and reliable contributions to the scientific community. [OMICS Group](#) hosts over 400 leading-edge peer reviewed Open Access Journals and organize over 300 International Conferences annually all over the world. OMICS Publishing Group journals have over 3 million readers and the fame and success of the same can be attributed to the strong editorial board which contains over 30000 eminent personalities that ensure a rapid, quality and quick review process.

About Omics Group conferences

- [OMICS Group](#) signed an agreement with more than 1000 International Societies to make healthcare information Open Access. [OMICS Group](#) Conferences make the perfect platform for global networking as it brings together renowned speakers and scientists across the globe to a most exciting and memorable scientific event filled with much enlightening interactive sessions, world class exhibitions and poster presentations
- Omics group has organised 500 conferences, workshops and national symposium across the major cities including San Francisco, Omaha, Orlando, Raleigh, Santa Clara, Chicago, Philadelphia, United Kingdom, Baltimore, San Antonio, Dubai, Hyderabad, Bangaluru and Mumbai.



OUTLINE

Laser Optics
2014

Philadelphia,
PA

September,
8-10, 2014

Mid-infrared semiconductor laser based trace gas analyzers: advances, applications & future outlook

S. So^{1a}, R. Lewicki^{1,2},
W. Ren², W. Jiang², Y. Cao², D. Jiang², F.K. Tittel^{2b}
V. Spagnolo³, Pietro Patimisco³

¹Sentinel Photonics, Princeton, NJ 08540

²Rice University, 6100 Main St., Houston, TX 77005

³University of Bari, Italy


a) sso@sentinelphotonics.com; b) fkt@rice.edu

- New Laser Based Trace Gas Sensor Technology
 - Novel Multipass Absorption Cell & Electronics
 - Quartz Enhanced Photoacoustic Spectroscopy
- Examples of Mid-Infrared Sensor Architectures
 - C₂H₆, NO, CO and CH₄
 - Future Directions of Laser Based Gas Sensor Technology and Conclusions

Research support by NSF ERC MIRTHE, NSF-ANR NexCILAS, the Robert Welch Foundation, and Sentinel Photonics Inc. via an EPA and NSF SBIR sub-award is acknowledged

Motivation for Mid-infrared C₂H₆ Detection

- Atmospheric chemistry and climate
 - Fossil fuel and biofuel consumption,
 - biomass burning,
 - vegetation/soil,
 - natural gas loss
- Oil and gas prospecting
- Application in medical breath analysis (a non-invasive method to identify and monitor different diseases):
 - asthma,
 - schizophrenia,
 - Lung cancer,
 - vitamin E deficiency.




HITRAN absorption spectra of C₂H₆, CH₄, and H₂O

NOAA Monitoring & Sampling Location: Alert, Nunavut, Canada



ALT, Ethane Concentration Measurements



General View on the Facility


Latitude: 82.4508° North

Longitude: 62.5056° West


Elevation: 200.00 m




C₂H₆ Detection with a 3.36 μm CW DFB LD using a Novel Compact Multipass Absorption Cell and Control Electronics



Schematic of a C₂H₆ gas sensor using a Nanoplus 3.36 μm DFB laser diode as an excitation source. M – mirror, CL – collimating lens, DM – dichroic mirror, MC – multipass cell, L – lens, SCB – sensor control board.



Innovative long path, small volume multipass gas cell: 57.6m with 459 passes




2f WMS signal for a C₂H₆ line at 2976.8 cm⁻¹ at a pressure of 200 Torr

Minimum detectable C₂H₆ concentration is:
~ 740 pptv (1σ; 1 s time resolution)




MPC dimensions: 17 x 6.5 x 5.5 (cm)
 Distance between the MPC mirrors: 12.5 cm


MULTIPASS CELL TECHNOLOGY




- High pathlength/volume ratio
- Simple spherical mirrors
- Utilize entire mirror surface - embrace optical aberration
- Flexible design – two opposing mirrors
- Typ. 10% throughput for 459 passes



Herriott Cell Pattern




Sentinel 3.7 m cell



Sentinel 57 m cell

From Conventional PAS to QEPAS



Quartz Tuning Fork as a Resonant Microphone for QEPAS




Unique properties

- Extremely low internal losses:
 - $Q \sim 10\,000$ at 1 atm
 - $Q \sim 100\,000$ in vacuum
- Acoustic quadrupole geometry
 - Low sensitivity to external sound
- Large dynamic range ($\sim 10^6$) – linear from thermal noise to breakdown deformation
 - 300K noise: $x \sim 10^{-11}$ cm
 - Breakdown: $x \sim 10^{-2}$ cm
- Wide temperature range: from 1.6K to ~ 700 K

Acoustic Micro-resonator (mR) tubes


- Optimum inner diameter: 0.6 mm; mR-QTF gap is 25-50 μm
- Optimum mR tubes must be ~ 4.4 mm long ($\sim \lambda/4 < l < \lambda/2$ for sound at 32.8 kHz)
- SNR of QTF with mR tubes: $\times 30$ (depending on gas composition and pressure)



Motivation for Nitric Oxide Detection



- Atmospheric Chemistry
- Environmental pollutant gas monitoring
 - NO_x monitoring from automobile exhaust and power plant emissions
 - Precursor of smog and acid rain
- Industrial process control
 - Formation of oxynitride gates in CMOS Devices
- NO in medicine and biology
 - Important signaling molecule in physiological processes in humans and mammals (1998 Nobel Prize in Physiology/Medicine)
 - Treatment of asthma, COPD, acute lung rejection
- Photofragmentation of nitro-based explosives

Molecular Absorption Spectra within two Mid-IR Atmospheric Windows and NO absorption @ 5.26 μm






Source: HITRAN 2000 database

Emission spectra of a 1900cm^{-1} TEC CW DFB QCL and HITRAN Simulated spectra



Output power: 117 mW @ 25 C
Thorlabs/Maxion

CW TEC DFB QCL based QEPAS NO Gas Sensor



CW HHL TEC DFB-QCL package and IR camera image of the laser beam at 630 mA and 20.5 deg C through tubes after ADM



Compact Prototype NO Sensor (September 2012)


Schematic of a DFB-QCL based Gas Sensor.

PcL – plano-convex lens, Ph – pinhole,


QTF – quartz tuning fork, mR – microresonator,

RC- reference cell, P-elec D – pyro electric detector

Performance of CW DFB-QCL based WMS QEPAS NO Sensor Platform



2f QEPAS signal (navy) and reference 3f signal (red) when DFB-QCL was tuned across **1900.08 cm⁻¹** NO line.





2f QEPAS signal amplitude for 95 ppb NO when DFB-QCL was locked to the **1900.08 cm⁻¹** line.

Minimum detectable NO concentration is:
~ 3 ppbv (1 σ ; 1 s time resolution)


Motivation for Carbon Monoxide Detection

- Atmospheric Chemistry
 - Incomplete combustion of natural gas, fossil fuel and other carbon containing fuels.
 - Impact on atmospheric chemistry through its reaction with hydroxyl (OH) for troposphere ozone formation and changing the level of greenhouse gases (e.g. CH₄).
- CO in medicine and biology
 - Hypertension, neurodegenerations, heart failure and inflammation have been linked to abnormality in CO metabolism and function.

Performance of a NWU 4.61 μm high power CW TEC DFB QCL




CW DFB-QCL optical power and current tuning at a four different QCL temperatures.




Estimated max wall-plug efficiency (WPE) is $\sim 7\%$ at 1.25A QCL drive-current.


CW DFB-QCL based CO QEPAS Sensor Results



2f QEPAS signal for dry (red) and moistured (blue) 5 ppm CO:N₂ mixture near 2169.2 cm⁻¹.



Dilution of a 5 ppm CO reference gas mixture when the CW DFB-QCL is locked to the 2169.2 cm⁻¹ R6 CO line.





Atmospheric CO concentration levels on Rice University campus, Houston, TX

Minimum detectable CO concentration is:
~ 2 ppbv (1 σ ; 1 s time resolution)

QEPAS based CH₄ and N₂O Gas Sensor

Motivation for CH₄ and N₂O Detection

- Prominent greenhouse gases
- Sources : Wetlands, leakage from natural gas systems, fossil fuel production and agriculture
- Applications: Environmental, medical and aerospace (N₂O)







Detection Limit (1 σ) with a **1-sec** averaging time

Methane (CH₄) (1275.04 cm⁻¹) **13 ppbv**

Nitrous Oxide (N₂O) (1275.5 cm⁻¹) **6 ppbv**

Deduced N₂O concentration in the ambient laboratory air: **331 ppbv**

QEPAS based CH₄ and N₂O Gas Sensor



QEPAS Sensor Control Board

QCL Current and TEC Driver, Performing wavelength modulation, Data acquisition, Applying continuous saw-tooth current ramping at 8 Hz, Testing QTF and low noise pre amplifier

QEPAS Performance for Trace Gas Species (September 2014)

	Molecule (Host)	Frequency, cm ⁻¹	Pressure, Torr	NNEA, cm ⁻¹ W/Hz ^{1/2}	Power, mW	NEC (τ=1s), ppmv
VIS	O ₃ (air)	35087.70	700	3.0×10 ⁻⁸	0.8	1.27
	O ₂ (N ₂)	13099.30	158	4.74×10 ⁻⁷	1228	13
	C ₂ H ₂ (N ₂)*	6523.88	720	4.1×10 ⁻⁹	57	0.03
NIR	NH ₃ (N ₂)*	6528.76	575	3.1×10 ⁻⁹	60	0.06
	C ₂ H ₄ (N ₂)*	6177.07	715	5.4×10 ⁻⁹	15	1.7
	CH ₄ (N ₂ +1.2% H ₂ O)*	6057.09	760	3.7×10 ⁻⁹	16	0.24
	N ₂ H ₄	6470.00	700	4.1×10 ⁻⁹	16	1
	H ₂ S (N ₂)*	6357.63	780	5.6×10 ⁻⁹	45	5
	HCl (N ₂ dry)	5739.26	760	5.2×10 ⁻⁸	15	0.7
	CO ₂ (N ₂ +1.5% H ₂ O) *	4991.26	50	1.4×10 ⁻⁸	4.4	18
	CH ₂ O (N ₂ :75% RH)*	2804.90	75	8.7×10 ⁻⁹	7.2	0.12
	CO (N ₂ +2.2% H ₂ O)	2176.28	100	1.4×10 ⁻⁷	71	0.002
	CO (propylene)	2196.66	50	7.4×10 ⁻⁸	6.5	0.14
Mid-IR	N ₂ O (air+5% SF ₆)	2195.63	50	1.5×10 ⁻⁸	19	0.007
	C ₂ H ₅ OH (N ₂)**	1934.2	770	2.2×10 ⁻⁷	10	90
	NO (N ₂ +H ₂ O)	1900.07	250	7.5×10 ⁻⁹	100	0.003
	C ₂ HF ₅ (N ₂)***	1208.62	770	7.8×10 ⁻⁹	6.6	0.009
	NH ₃ (N ₂)*	1046.39	110	1.6×10 ⁻⁸	20	0.006
	SF ₆	948.62	75	2.7×10 ⁻¹⁰	18	5×10 ⁻⁵ (50 ppt)

* - Improved microresonator

** - Improved microresonator and double optical pass through ADM



*** - With amplitude modulation and metal microresonator

NNEA – normalized noise equivalent absorption coefficient.


NEC – noise equivalent concentration for available laser power and τ=1s time constant, 18 dB/oct filter slope.

For comparison: conventional PAS $2.2 \times 10^{-9} \text{ cm}^{-1}\text{W}/\sqrt{\text{Hz}}$ for NH₃

Mini Methane Sensor for UAVs





Miniaturization



Remote controlled quad-copter UAV
for pipeline sniffing -
payload maximum only 600g!

Mini Methane Sensor for UAVs

Sensor Performance



Onboard pressure controller and pump system
Continuous measurement [10Hz] with onboard processing
Direct concentration output – no post-processing necessary



Future Directions and Outlook

- New target analytes such as carbonyl sulfide (OCS), formaldehyde (CH₂O), nitrous acid (HNO₂), hydrogen peroxide (H₂O₂), ethylene (C₂H₄), ozone (O₃), nitrate (NO₃), propane (C₃H₈), and benzene (C₆H₆)
- Ultra-compact, low cost, robust sensors (e.g. C₂H₆, NO, CO...)
- Monitoring of broadband absorbers: acetone (C₃H₆O), acetone peroxide (TATP), UF₆...
- Optical power build-up cavity designs
- Development of trace gas sensor networks
- QEPAS based detection at THz frequencies

Future Directions and Outlook

- Development of robust, compact sensitive, mid-infrared trace gas sensor technology based on room temperature, continuous wave, DFB QCL and ICLs for environmental, industrial, biomedical monitoring and security applications
- Seven target trace gas species were detected with a 1 sec sampling time:
 - C_2H_6 at $\sim 3.36 \mu\text{m}$ with a detection sensitivity of 740 pptv using TDLAS
 - NH_3 at $\sim 10.4 \mu\text{m}$ with a detection sensitivity of ~ 1 ppbv (200 sec averaging time);
 - NO at $\sim 5.26 \mu\text{m}$ with a detection limit of 3 ppbv
 - CO at $\sim 4.61 \mu\text{m}$ with minimum detection limit of 2 ppbv
 - SO_2 at $\sim 7.24 \mu\text{m}$ with a detection limit of 100 ppbv
 - CH_4 and N_2O at $\sim 7.28 \mu\text{m}$ currently in progress with detection limits of 20 and 7 ppbv, respectively.
- New target analytes such as CH_2O , H_2O_2 , and C_2H_4 ,
- Monitoring of broadband absorbers such as acetone, C_3H_8 , C_6H_6 and UF_6



Let Us Meet Again

We welcome all to our future group conferences of Omics group international

Please visit:

www.omicsgroup.com

www.Conferenceseries.com

<http://optics.conferenceseries.com/>